Product Specification Microlite[®]





Trade name: Microlite[®] Intended use: To protect the user from harmful chemical and substance contamination within food and medical industry and trade based work environments. Materials: 100% synthetic nitrile Colour: Blue. Type: Nitrile Examination Glove, Powder Free, Online Single Chlorinated, Non-sterile, Ambidextrous. Surface finish: Fully textured, beaded cuff. Powder: No powder lubricant added. Storage condition: The gloves shall maintain their properties when stored in a dry condition. Avoid direct sunlight. Shelf life: The gloves shall have shelf life of 5 years from the date of manufacture with the above storage conditions. UOM: 300 gloves. Small, Medium, Large, XL. Sizes:

Medical Standards

Microlite® meets ATSM D6319 for Medical Exam use



Product Specification

Dimensions	Size	Palm width	Length (mm)
Millimetres (Inches)	Small	85 ± 5mm (3.34" ± 0.1")	240 (minimum) (9.44")
	Medium	95 ± 5mm (3.74" ± 0.1")	240 (minimum) (9.44")
	Large	106 ± 5mm (4.17" ± 0.1")	240 (minimum) (9.44")
	X Large	113 ± 5mm (4.45" ± 0.1")	240 (minimum) (9.44")

Thickness	Palm	Fi
Millimetres	>0.05mm	>(
(Inches)	(>0.00197")	(>(

inger 0.06mm 0.00236")

Box



Quality specification* And acceptable

number of defects.

Inspection criteria Dimension

F

Physical properties Water tight test Visual inspection - major Visual inspection - minor **Residual powder** Mix size/ mix glove/ mix ha

Packaging Configuration & Labelling:

Packaging style:	300 gloves X 10 boxes X 1 carton
Size marking:	The size of the gloves shall be marked in the check box on every carton with black ink.
Packing:	300 gloves per inner box (by weight) 10 boxes per outer carton.
Labelling:	The inner boxes and outer cartons have the following information: product description, brand, size and net quantity of gloves, country of origin, name and address of distributor, and single use only symbol.
Barcodes:	GS1 format.

	Sampling Plan/ AQL
	S2 - AQL 4.0
	S2 - AQL 4.0
	G1 - AQL 1.5
	G1 - AQL 2.5
	G1 - AQL 4.0
	n = 5 (average)
and	Not allowed
	l de la constante de



For more information: customerservice@theglovecompany.com

1st Oct 2018